



GUY COTTEN

150 N LIFE JACKET PERF



Technical data

PROTECTION FRAMEWORK

ISO 12402-3: 2006 – Life jackets, performance level 150

A life jacket that will lift an exhausted or unconscious man overboard to the surface as quickly as possible, without affecting his health and turn him in a position that will allow him to breath while waiting to be rescued.

ISO 12401: 2009 – Deck safety harness and safety lines on recreational craft

Intended to be worn by all persons exposed by their work to the risk of falling overboard and drowning.

Universal use: for fishing, coastal waters and open seas.

PERFORMANCE



- Buoyancy of around 170 Newtons
- Completely automatic self-inflation
- Manual trigger device included
- CO2 indicator according to the depth of thread
- Content of the gas cartridge = 32 g CO2
- Quick and easy pill replacement, cartridge loading and rearming
- The buoyancy chamber is covered with a 300-cm² retroreflective strip
- Deck harness included

Outer envelope made of Nylpêche fabric, withstanding wear better



It is worn high on the body, which prevents wear, and movement is not impaired.

Easy to put on with a stainless-steel buckle fastener at the front.

Lung with double wall

Manual pull cord trigger and mouthpiece for oral inflation

2 tone whistle

Equipped with a fixation for a multi-flash light or beacon



Option: MANUAL KIT overriding self-inflation

Option: Can be worn under the 3D vest or under the ISOTOP smock

PRINCIPAL MATERIALS

- Polyamide
- Polyurethane coating

i For maximum protection, please check the condition of your life jacket, the trigger devices and its servicing date.

To protect the trigger bobbin/pill from dissolving prematurely due to high temperatures/atmospheric humidity, adapt the regularity with which you change these bobbins/pills.

Annual servicing is a must. Have your jacket serviced by one of the approved GUY COTTEN servicing stations.

A neoprene cuff (sold on option) can be adapted to prevent the buckle catching.

INSTRUCTIONS

- CARE**
 - Regularly use oral inflation to check the life jacket.
- STORAGE**
 - 54°C to +70°C
 - Never fold it when wet.
 - Before storing, check the jacket is properly rearmed (cf Kit)

Colour: yellow

Tailles: one size

Weight of person: 40 to 120 kg

Weight of jacket: about 1.5 kg

The lifespan of your equipment is linked to conditions of use and quality of care